





10X

- History of STx
- Current STx technologies
 - 10X Genomics Visium
 - NanoString GeoMx


RNA SEQUENCING IN SITU











ANALYTICAL CHALLENGES OF SPATIALLY-RESOLVED 'OMICS DATA









10X

GMX



How to leverage the spatial aspect of the data effectively?



SPATIAL DATA ANALYSIS IN MOLECULAR BIOLOGY



Data with a spatial dimension have distinct properties




A different set of statistical and inferential considerations.



10X

GMX



Problem

Regression

Coefficient


APPLICATION OF GEOSPATIAL METHODS IN BIOLOGY





10X

Space has the **same features** in all disciplines



SPATIAL TRANSCRIPTOMICS IS NOT NEW... 1989 WMISH in Drosophila 1969 Radioactive ISH of rRNA 1991 in situ 1977 FISH of rRNA 1982 FISH of actin reporter in C. mRNA elegans 1969 1971 1975 1979 1981 1985 1987 1991 1977 983 1989 1982 Immunological 1987 1973 Radioactive FISH with biotin Drosophila ISH of globin mRNAs labeled probe enhancer trap 1989 ES cell enhancer and gene trap in mice 2013 High Prequel Type throughput RCA + ISS 1996 Commercial LCM 2012 Tomo-array of 1988 Ligase SNV 1995 cDNA 2002 Combinatorial mouse brain FISH for mRNA detection microarray 2015 MERFISH 1976 1978 1980 1982 1984 1986 1988 1990 1992 1994 1996 1998 2000 2002 2004 2006 2008 2010 2012 2014 2016 2018 2020 1998 Solexa founded 2019 GeoMX DSP 1989 Single cell 1976 LCM 2008 RNA-seq cDNA amplification 1999 LCM + 2014 seqFISH microarray 1989 FISH with 1998 smFISH combinatorial 2016 Spatial Transcriptomics barcoding

Type

© ICBAM research group, Newcastle University, UK

ROI selection

NGS barcoding

smFISH

Source: Museum of spatial transcriptomics (https://doi.org/10.1038/s41592-022-01409-2)

10X

GMX

_/

CURRENT METHODS OF SPATIAL TRANSCRIPTOMICS





10X

GMX

_			
Technology	Highest	Max Number	Number of unique
	Resolution	of Targets	spatial barcodes
smFISH	single cell	10	N/a (probe-based)
seqFISH	single cell	10421	N/a (probe-based)
seqFISH+	single cell	10000	N/a (probe-based)
MERFISH	single cell	4209	N/a (probe-based)
STARmap	single cell	1020	N/a (probe-based)
MOSAICA	single cell	10	N/a (probe-based)
original ST	100μm	transcriptome	1007
10X Visium	55μm	transcriptome	~5000
HDST	2μm	transcriptome	~1.54 million
Slide-seq	10μm	transcriptome	10000
Slide-seqV2	10μm	transcriptome	10000
DBiT-seq	10μm	transcriptome	2500
PIXEL-seq	1.22μm	transcriptome	N/a (Illumina flow-like)
space-TREX	55μm	transcriptome	~5000
sci-Space	73.2μm	transcriptome	7506
GeoMx DSP	5μm	20175	N/a (probe-based)
Stereo-seq	500nm	transcriptome	Approx. 800m


Current spatial transcriptomics technologies and their main specifications.

ARRAY AND PROBE EXAMPLES










Array examples

10X

GMX



DNA oligo tag UΥ photocleavable linker **DNA** probe **RNA**

Probe examples

10X Visium

GeoMx




Visium Slide




Sample Prep & Staining

10X






10X




GMX

Fluorescent 10X




H&E





Visium Slide



Sample Prep & Staining


Imaging

Permeabilisation & RNA Diffusion

10X




Visium Slide



10X



Visium Slide




10X





10X




GMX

Tissue on slide





10X



GMX


Tissue on slide

Probes incubation






10X








10X





Eleftherios Zormpas





Dr Simon J Cockell






Dr Rachel Queen





Prof. Alex Comber





iSMB feedback form:









